

## **PRODUCT DATA SHEET**

## Anti-Flotillin-1 Rabbit Polyclonal antibody

Cat Code: HBM-FLOT-SR1-###

Size: 100 or 50 µg

Description	Flotillin belongs to the band 7.2/stomatin protein family and appears to be strongly expressed in muscle cells and fibroblasts. Flotillin expression is also correlated with Alzheimer development. Flotillin is highly expressed on exosomes and appears to be involved in exosome release
MW	Flotillin-1MW: 47 kDa
Immunogen	aa 324-427 of Flotillin-1 of human origin
Туре	Rabbit Polyclonal
Species Reactivity	Human.
Concentration	0,2 mg/ml.
Application	ELISA, Flow Cytometry, Western blotting
Recommended dilution	Recommended dilution for ELISA and WB 1/500
Storage buffer	Phosphate buffered saline (PBS) with 15 mM sodium azide, pH 7.4
Storage Condition	Store at 2°-8°C. Do not freeze. Do not use after expiration date indicated on vial label
Expiration	See vial label
Lot Number	See vial label
Product notes	Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing

This product is for research use only and is not approved for use in humans or in clinical diagnosis.

HansaBioMed LLC Akadeemia tee 15, 12618 Tallinn, ESTONIA

## www.hansabiomed.eu Tel: +372 6204348 Fax: +372 6204358 email: info@hansabiomed.eu



Contact in USA Galen Laboratory Supplies P.O. Box 682 Middletown, CT 06457. USA

www.galenlabsupplies.com Tel: 860-704-9058 Fax: 888-260-6091 email: info@ galenlabsupplies.com